

Summary of the Spectra Experiment

This document summarises the *Spectra Experiment*. The full description of the experiment appears at:

<http://www.upscale.utoronto.ca/IYearLab/Intros/Spectra/Spectra.html>.

Apparatus Notes

- Align the cross hairs in the eyepiece to form an x, not a +.
- Move the eyepiece in and out until the focus of the cross hairs is as sharp as possible.
- Adjust the slits to just resolve the two yellow lines in the Sodium Spectrum.

Calibration

Use the Helium lamp for calibration. Fit your data to the Hartmann dispersion relation:

$$y = \frac{m}{\lambda - \lambda_0} + b$$

λ_0 is a constant for the spectrometer, and is supplied.

Hydrogen Spectrum

Measure the wavelengths of the lines in the Hydrogen spectrum to determine R_H where:

$$\frac{1}{\lambda} = R_H \left(\frac{1}{2^2} - \frac{1}{n^2} \right), n = 3, 4, 5, \dots$$

Significance of R_H

Calculate $h c R_H$, where h is *Planck's constant*, and c is the speed of light. Convert to *electron-volts* and compare to the ionization energy of atomic Hydrogen, which is 13.6 eV.

Use the equation in *Preparatory Question 4* to discuss the significance of your result.

Gas Identification Sleuthing

Choose at least one of the unknowns and identify it.

Preparatory Questions

- For the figures of the Scale and the Vernier in the full description of the experiment, what is the reading error? Note that the answer will depend on at least:
 - Your vision.
 - The quality of the computer monitor you are using.
 - The resolution of the figures themselves (90 x 90 pixels).

For the real apparatus, the reading error will depend on at least your vision and the quality of the engraving of the lines onto the spectrometer.

- You measure the position of a spectral line, and get a scale reading of 9.26. You decide that the reading error is 0.05. Thus:

$$y = 9.26 \pm 0.05$$

Your calibration of the spectrometer using Equation 1 resulted in the values:

$$b = 4.1275 \pm 0.0032$$

$$m = 1491.2 \pm 3.5 \text{ nm}$$

Your spectrometer has a value $\lambda_0 = 283.2 \pm 0.4 \text{ nm}$.

- What is the value and error of $(y - b)$? What is the dominant error in this value?
- What is the value and error of $m/(y - b)$? What is the dominant error in this value?
- What is the value and error of $m/(y - b) + \lambda_0$? What is the dominant error in this value? Note that this is the value of the wavelength.

Note that all the above numbers are fictitious. However, the principles of error propagation and especially learning to ignore non-dominant errors will be the same for your real data.

- When you measure the Hydrogen spectrum, you will see only four or possibly five lines. However, the Balmer formula seems to predict an infinite number of lines corresponding to the infinite number of integers greater than or equal to 3. Do you think those lines exist in the spectrum?
 - If your answer is *Yes* why don't you see them?
 - If your answer is *No* why don't they exist?

- A simple variation of the Balmer equation is:

$$\frac{1}{\lambda} = R_H \left(\frac{1}{1^2} - \frac{1}{n^2} \right), n = 2, 4, 5, \dots$$

Do you think these lines exist in the Hydrogen spectrum?


- If your answer is *Yes* why don't you see them?
- If your answer is *No* why don't they exist?

Spectral Wavelength Tables

- Use the intensity indications with caution. They are only a general guide, and your lines may have different intensities.

- The tables give most of the lines you will be able to see, and many that you won't be able to see if you are using a narrow slit width. However, they are not complete.
- Lines separated by less than 1 nm will not be resolved if the slit is wide. If the slit is narrow, weak lines won't be seen.
- You may assume that errors in the wavelengths are negligible. Typically wavelengths are known to 0.00001 nm or better.
- A full set of wavelength tables is maintained by the U.S. National Institute for Standards and Technology at http://physics.nist.gov/cgi-bin/AtData/lines_form.

<u>HELIUM</u>		
<i>WAVELENGTH</i> <i>nm</i>	<i>RELATIVE</i> <i>INTENSITY</i>	<i>COLOUR</i>
728.1	2	RED
706.5	4	RED
667.8	6	RED
656.0	1	RED
587.6	10	YELLOW
504.8	4	GREEN
501.6	6	GREEN
492.2	5	GREEN
485.9	2	GREEN
471.3	5	BLUE
447.1	6	BLUE
443.8	1	VIOLET
438.8	4	VIOLET
416.9	1	VIOLET
414.4	2	VIOLET
412.1	3	VIOLET
402.6	4	VIOLET
396.5	1	VIOLET
388.9	3	VIOLET



<u>ARGON</u>			<u>KRYPTON</u>		
<i>WAVELENGTH</i> nm COLOUR	<i>RELATIVE</i> INTENSITY	COLOUR	<i>WAVELENGTH</i> nm	<i>RELATIVE</i> INTENSITY	
574.0	2	GREEN	645.6	5	RED
565.0	3	GREEN	642.1	5	RED
560.7	3	GREEN	605.6	2	RED
557.3	3	GREEN	601.2	2	ORANGE
549.6	3	GREEN	599.4	2	ORANGE
522.1	2	GREEN	588.0	1	ORANGE
518.8	3	GREEN	587.1	10	ORANGE
516.2	3	GREEN	584.1	1	YELLOW
470.2	1	BLUE	583.3	1	YELLOW
462.8	1	BLUE	570.8	1	GREEN
459.6	1	BLUE	567.2	1	GREEN
452.2	1	VIOLET	565.0	1	GREEN
451.1	2	VIOLET	558.0	1	GREEN
433.5	2	VIOLET	557.0	10	GREEN
433.4	2	VIOLET	556.2	2	GREEN
430.0	3	VIOLET	450.2	5	VIOLET
426.6	3	VIOLET	446.4	5	VIOLET
425.9	3	VIOLET	445.4	5	VIOLET
420.1	2	VIOLET	440.0	2	VIOLET
419.8	2	VIOLET	437.6	5	VIOLET
416.4	3	VIOLET	436.3	4	VIOLET
415.9	2	VIOLET	432.0	3	VIOLET
431.9	2	VIOLET	427.4	5	VIOLET

Argon has many faint lines in the red and yellow which vary in intensity depending on the source and because of the confusion that this can lead to only wavelengths less than 580 nm are given. In this region there are a very large number of lines. Only relatively brighter ones are listed. Fainter ones may provide a haze in the background.

<u>MERCURY</u>			<u>XENON</u>		
<i>WAVELENGTH</i>	<i>RELATIVE</i>	<i>COLOUR</i>	<i>WAVELENGTH</i>	<i>RELATIVE</i>	<i>COLOUR</i>
<i>nm</i>	<i>INTENSITY</i>		<i>nm</i>	<i>INTENSITY</i>	
708.2	1	RED	647.3	2	RED
704.5	2	RED	647.0	3	RED
690.7	1	RED	631.8	5	RED
671.6	1	RED	620.1	1	RED
658.5	1	RED	619.8	1	RED
638.3	2	RED	618.2	3	RED
623.4	2	RED	618.0	1	RED
612.3	2	RED	617.8	2	RED
607.3	2	ORANGE	616.4	1	RED
602.4	2	ORANGE	593.4	2	ORANGE
601.7	1	ORANGE	593.1	1	ORANGE
589.0	1	YELLOW	589.5	2	ORANGE
579.1	8	YELLOW	587.5	1	ORANGE
577.0	6	YELLOW	582.5	2	YELLOW
567.7	1	YELLOW	582.4	3	YELLOW
567.6	1	YELLOW	571.6	1	YELLOW
546.1	10	GREEN	569.7	1	YELLOW
536.5	1	GREEN	569.6	1	YELLOW
520.5	1	GREEN	546.0	1	GREEN
519.6	1	GREEN	539.3	1	GREEN
512.1	1	GREEN	502.8	3	GREEN
504.6	1	GREEN	492.3	4	GREEN
502.6	1	GREEN	491.7	4	GREEN
496.0	1	GREEN	484.3	4	GREEN
491.6	5	BLUE	483.0	4	GREEN
452.3	1	BLUE	480.7	5	GREEN
435.8	6	VIOLET	479.3	1	BLUE
434.8	2	VIOLET	473.4	5	BLUE
433.9	1	VIOLET	469.7	4	BLUE
421.2	1	VIOLET	467.1	10	BLUE
420.6	1	VIOLET	462.4	5	BLUE
415.7	1	VIOLET	458.3	1	VIOLET
407.8	5	VIOLET	452.5	2	VIOLET
414.7	5	VIOLET	450.1	2	VIOLET

<u>NEON</u>			<u>NITROGEN</u>		
<i>WAVELENGTH</i> <i>nm</i>	<i>RELATIVE</i> <i>INTENSITY</i>	<i>COLOUR</i>	<i>WAVELENGTH</i> <i>nm</i>	<i>RELATIVE</i> <i>INTENSITY</i>	<i>COLOUR</i>
724.5	1	RED	497.6	3	GREEN
717.4	1	RED	491.7	3	GREEN
703.2	5	RED	481.5	3	GREEN
702.4	3	RED	472.4	3	BLUE
692.9	6	RED	466.7	3	BLUE
667.8	7	RED	464.9	3	BLUE
659.9	7	RED	457.4	3	BLUE
653.3	7	RED	449.0	3	BLUE
650.7	7	RED	441.7	3	VIOLET
609.6	5	ORANGE	435.5	3	VIOLET
607.4	7	ORANGE	434.4	2	VIOLET
603.0	5	ORANGE	427.0	2	VIOLET
596.5	4	ORANGE	420.1	2	VIOLET
588.2	6	YELLOW	414.2	2	VIOLET
585.2	10	YELLOW	409.5	1	VIOLET
540.1	5	GREEN	406.0	1	VIOLET
			399.8	1	VIOLET
			394.3	1	VIOLET

Many orange and yellow lines have been omitted as well as all lines of wavelength less than 540 nm (hundreds). Most of these are faint but some overlap gives the appearance of bright lines.

Since **nitrogen** is a molecule, the spectrum consists of bands rather than lines. This is due to rotation of the molecules. In the visible the most prominent structure is the First Positive series with about 30 regular spaced bands in the region 500-700 nm. Only the band heads of the Second Positive series are tabled above. The bands trail off to shorter wavelengths. As indicated by the relative intensities on a scale of 10, the Second Positive series is less intense than the First Positive series.